


INFUSION/BLOOD WARMER

ANIMEC AM-301

RELIABLE, SAFE, COST-EFFECTIVE

NO DISPOSABLES
ADJUSTABLE TEMPERATURE
TEMPERATURE DISPLAY
OVER-TEMPERATURE ALARM
SIMPLE OPERATION

With the ANIMEC AM-301, warmed fluid can be efficiently delivered to patients with ease, safety and no risk of contamination. The ANIMEC AM-301 will surely benefit your patients and hospital.



ELLTEC

Hypothermia, which may result in serious side-effects, has always been a major concern. Infusion/ Transfusion at body temperature helps prevent and treat hypothermia and also assists effective intravenous administration.

The **ANIMEC AM-301** has various special features to

- **meet individual needs**

User-adjustable set-point temperature
[36°C or 39°C]

- **ensure safe use**

Easy to read LED display of actual outlet fluid temperature
Audible and Visible over-temperature alarms
Alarm test available by the "check" button

- **save cost**


No disposable set required
Available for two standard tube sizes
(4mm & 5mm)

- **extend application**

User-selectable flow rate range [gravity drip]
1-15ml/min: Single-channel setting [photo1]
1-20ml/min: Dual-channel setting [photo2]

- **offer greater convenience**

Mountable horizontally and vertically
Compact and light-weight with a handle
Easy maintenance



"Simple operation/Quick warm up"

1. Set a tube into the channel and fill it up with fluid.
2. Turn on the power.
3. Adjust set-point temperature and push the "start" button.

APPLICATIONS

Infusion or Transfusion for pre-surgery,
post-surgery or short-term surgery

SPECIFICATIONS

| Type | AM-301-4A/5A | AM-301-4B/5B |
|-----------------------|------------------------------------|------------------------------------|
| Tube Size | 4A: 3.1-4.0mm φ 5A: 4.1-5.0mm φ | 4B: 3.1-4.0mm φ 5B: 4.1-5.0mm φ |
| Flow Rate | 1 - 20 ml/min | |
| Set-point Temperature | 36°C / 39°C | |
| Alarm Temperature | 42.5°C (Outlet Fluid Temperature) | |
| Heat Exchanger | Dry Heat [Silicone Rubber Heater] | |
| Power Requirements | 220-250V~ 50/60Hz | 110-120V~ 50/60Hz |
| Dimensions | L 115 x W 65 x H 265 (mm) | |
| Weight | 1 kg | |

For further information, please contact our distributor or Elltec on E-mail:international@elltec.co.jp

ELLTECco.,LTD.

Shirakawa No.6 Bldg. 2-18-5, Nishiki Naka-ku, Nagoya JAPAN
Tel: +81-52-201-7308 Fax: +81-52-232-2870
E-mail:international@elltec.co.jp